



Technical Data Sheet : IPA30 and IPA40



Technical Specifications	IPA 30	IPA 40
Type of Machine	Automatic Mounting Press with 30 mm dia. mold size	Automatic Mounting Press with 40 mm dia. mold size
System and control	Pneumatic cylinder Microprocessor, CE safety rules	Pneumatic cylinder Microprocessor, CE safety rules
Whole structure	Oven painted steel	Oven painted steel
Pressure control	Manometer Restored automatically, 8 bar	Manometer Restored automatically, 8 bar
Heating system	Electric immersion with digitally control	Electric immersion with digitally control
Water cooling system	Automatic end of heating process	Automatic end of heating process
Heating and cooling temperature	Preset / Program key	Preset / Program key
Audio signal	Cycle is completed	Cycle is completed
Max. heating temperature and time	200° C / 99 min	200° C / 99 min
Warning lights	Heating / Cooling Set point 1 / Set point 2	Heating / Cooling Set point 1 / Set point 2
Visual display unit	Set parameters Real temperature during process	Set parameters Real temperature during process
Power	600 W, 220V 1-phase	650 W, 220V 1-phase
Weight	40 Kg	42 Kg
Dimension (mm) : W	400	400
D	400	400
H	510	510

REMET S.A.S.

Via Scarlatti 2, 40033 Casalecchio di Reno (Bologna) ITALY