



Technical Data Sheet : LS250V and LS250V TWIN



Technical Specifications	LS250V	LS250V TWIN
Type of Machine	Grinder/Polisher	Grinder/Polisher
Disc type and diameter of disc	Single disc 250 mm dia.	Double disc 250 mm dia.
Rotational speed	Variable speed 0 to 300 rpm Digital display	Variable speed 0 to 300 Digital display (OPT) separate control for each disc
Whole structure	Oven painted acid-proof steel	Oven painted acid-proof steel
Feet	Shock-isolating	Shock-isolating
Transmission of working disc	Precise and noiseless gear-reducer	Precise and Noiseless gear-reducer
Drainer	Polished stainless steel	Polished stainless steel
Adjusting the number of motor revolution	Potentiometer	Potentiometer
Water faucet	1 pc.	2 pcs.
Water nozzle	1 pc.	2 pcs.
Splash ring	Option	Option
Motor power	250 W	2x250 W
Power supply	220V 1-phase	220V 1-phase
Weight	32 Kg	60 Kg
Dimension (mm) : W	370	720
D	500	500
H	300	300

REMET S.A.S.

Via Scarlatti 2, 40033 Casalecchio di Reno (Bologna) ITALY